

Making *E. coli* Competent cells: (Rubidium chloride method)

Procedure:

1. Inoculate 1 ml from overnight culture into 100 ml LB broth (scale up or down as needed). Incubate at 37°C with aeration to A600=0.4
2. Keep on ice 15 min.
3. Pellet cells in 50 ml falcon tube (3-5000 x g 5 min)
4. Discard supernatant and add 1/3 volume of TfbI, resuspend and incubate back on ice for 15 min.
5. Pellet cells at low speed. (3000 x g 5 min).
6. Discard supernatant and resuspend in 1/30 volume of TfbII (1/30 of the culture volume), Keep on ice for 15 min and either use immediately or quick freeze at -70C for storage. I usually save these in 0.1 to 0.2 ml aliquots. Quick freeze in liquid nitrogen prior to storage in a -70 C freezer. Thaw on ice just before using in a transformation experiment.

I typically transform 100 ul cells with 2-10 ul of a ligation reaction, and you should get between 1×10^8 to 1×10^9 cfu's/ug DNA.

TFB I (per 200 ml)		
compound	amount	final conc.
potassium acetate	.588 g	30 mM
rubidium chloride	2.42 g	100 mM
calcium chloride	0.294 g	10 mM
manganese chloride	2.0 g	50 mM
glycerol	30 ml	15% v/v
<i>pH 5.8 with dilute acetic acid</i>		

TFB II (per 100 ml)		
compound	amount	Final conc.
MOPS	0.21 g	10 mM
calcium chloride	1.1 g	75 mM
rubidium chloride	0.121 g	10 mM
glycerol	15 ml	15% v/v
<i>pH 6.5 with dilute NaOH</i>		